


PRELIMINARY, FINAL IMPRESSION & BITE REGISTRATION

Parkell’s proficiency in manufacturing VPS impression and bite materials extends beyond 6 decades. As the market leader, Blu-Mousse® is still the #1 preference with clinicians worldwide. Cinch®, Sharp™, and Mach™ materials were developed with the same expertise as Blu-Mousse®. That’s why Parkell continues to be a leader in developing and manufacturing impression materials.



Cinch® Cartridges

- Exceptionally accurate, hydrophilic and strong.
- Use as a final impression material or an alginate substitute.
- Choice of viscosities and setting times.
- 5 Minute Quick-Pour Formula.



Name	Quick Cinch® Light	Cinch® Light	Cinch® 90 Medium	Cinch® Platinum Medium	Quick Cinch® Heavy	Cinch® AutoMix Putty
Stock No.	S461S	S460S	S453S	S454S	S463S	S456S
Set time	90 sec.	4 min.	90 sec	4 min.	90 sec.	4 min.
Working time	60 sec.	3 min.	60 sec.	3 min.	60 sec.	3 min.

Includes: Four 50ml cartridges (200ml ttl.) and mixing tips.



Cinch® Platinum Tubes

- No dispensing gun required.

Includes:
Four tubes (2 Catalyst & 2 Base) 320ml ttl.

Stock No.	Set Time	Working Time
S435	4 min.	3 min.



Cinch® 90 Super Cartridge

- Ideal for bulk mixing.

Includes:
One 380ml cartridge and mixing tips.

Stock No.	Set Time	Working Time
S464	90 sec.	60 sec.